Using the Alamut[®] Visual Plus Extranet

Accessing the Extranet

• As an Alamut[®] Visual Plus administrator for your institution, you will receive an email with a link allowing you to validate your administrator account and to choose a password.

Your username is embedded in this email.

 After having confirmed your account, the Extranet will be accessible using your institution ID, username and password. The URL is https://extranet.interactive-biosoftware.com.

	alamut VISUAL PLUS™
-	Alamut [®] Visual Plus extranet
Please S	ign In
Please S	-
	n ID
Institutio	m ID
İnstitutio User Na	m ID
İnstitutio User Na	m ID

Users Management

• Select "Users Management" in the menu.



- To add a user click the "Add user" button Add user
- A pop-up window will appear, where you can fill-in the details of the new user. This will automatically send an email to the new user prompting them to create a password to access the extranet.

• The list of users appears in a table:

Username \$	Name 🌻	First Name 🌲	Dpt \$	Created \$	Admin \$	Active \$	Connected \$	
			info	2017-12-08 09:29:39	~	~		c i
			Sales	2017-02-17 09:28:30		~		C i
				2018-01-24 09:27:49	~	~		C i
			Bioinfo	2018-02-05 09:50:33		~		C i
			Info		~	~		c i
			Sales	2017-12-27 16:53:24		-		c i
			Bioinfo	2018-01-22 10:03:39	~	×		c i
			Info			~		c i
			Sales	2018-01-24 12:01:09		~		c i
			П	2017-02-07 16:53:42	~	~		G

- You can edit a user's information by clicking on its "Edit User" button
- You can delete a user by clicking on its "Delete User" button 🔳 (this action is irreversible).

Activity statistics

• Select "Activity Stats" in the menu.



• In the "Parameters" section, click on the "Activity period" field to choose the time interval you want statistics on, and click on the "Apply" button.

	with t	perio	d:											
Ê	Octo	ober	19, 20	018 -	Octo	ber 22	, 2018	•						
1	0/19/	2018					1	0/22/	2018	3				Today
<		0	ct 20 ⁻	18					N	ov 20	18			Yesterday
Su	Мо	Tu	We	Th	Fr	Sa	Su	Мо	Tu	We	Th	Fr	Sa	Last 7 Days
30	1	2	3	4	5	6	28	29	30	31	1	2	3	Last 30 Days
7	8	9	10	11	12	13	4	5	6	7	8	9	10	This Month
14	15	16	17	18	19	20	11	12	13	14	15	16	17	Last Month
21	22	23	24	25	26	27	18	19	20	21	22	23	24	Custom Range
28	29	30	31	1	2	3	25	26	27	28	29	30	1	Apply Cancel
4	5	6	7	8	9	10	2	3	4	5	6	7	8	Apply

• The "Global activity" section shows the consolidated use of your Alamut[®] Visual Plus license for the chosen period.

Global activity						
		Global	activity of license	1000		
20						
10						
10						
0 Thu 00:00 Max tokens	Thu 12:00	Fri 00:00	Fri 12:00	Sat 00:00	Sat 12:00	Sun 00:0

• The "Activity by user" section shows the user-by-user consumption of your Alamut[®] Visual Plus license for the chosen period.

			Planar			
0 Thu 00:00	Thu 12:00	Fri 00:00	Fri 12:00	Sat 00:00	Sat 12:00	Sun 00:00
			Genes			
1 0 Thu 00:00	Thu 12:00	Fri 00:00	Fri 12:00	Sat 00:00	Sat 12:00	Sun 00:00
			Reality			
0 Thu 00:00	Thu 12:00	Fri 00:00	Fri 12:00	Sat 00:00	Sat 12:00	Sun 00:00
			Daniel			
0 Thu 00:00	Thu 12:00	Fri 00:00	Fri 12:00	Sat 00:00	Sat 12:00	Sun 00:00

License details

• Select "License Details" in the menu.



• In this section you can find all the necessary information about your Alamut[®] Visual Plus license.

cense details		
Institution ID	ATR1-0	
License key	11/2002/200	
Creation date	2018-07-02	
Expiration date	2019-07-02	
Token(s)	20	

• This also gives the administrator the ability to set a timeout option in the "Timeout Settings". This enables Alamut[®] Visual Plus to timeout after a set period of inactivity, say 15 minutes, freeing the license token for another user, but also saving the data of the inactive user so they can return to their analysis without disruption.

O Timeout settings	
Activate timeout	
Timeout starts after (minutes)	25
Timeout period lasts for (minutes)	30
	Save